



Mako Total Hip array instrument kit



**Upper tray**


| Description |                               | Description |                               |
|-------------|-------------------------------|-------------|-------------------------------|
| 1           | 112240 Pelvic array adapter   | 5           | 206693 MICS adapter cover     |
| 2           | 112390 3-Pin clamp            | 6           | 160230 RIO registration array |
| 3           | 112270 2-Pin clamp*           | 7           | 160245 Hip femoral array post |
| 4           | 204980 MICS reamer attachment | 8           | 160336 Hip probe              |

\*Optional, available upon request



**Lower tray**

| <b>Description</b> |                                   | <b>Description</b> |                                     |
|--------------------|-----------------------------------|--------------------|-------------------------------------|
| <b>1</b>           | 112227 Rio base array             | <b>5</b>           | 110755 Postero lateral broach array |
| <b>2</b>           | 206967 or 205020 Hip end effector | <b>6</b>           | 110765 Antero lateral broach array  |
| <b>3</b>           | 112230 Pelvic array               |                    |                                     |
| <b>4</b>           | 112330 Hip femoral array          |                    |                                     |




**Description**

|           |             |                                   |
|-----------|-------------|-----------------------------------|
| <b>1</b>  | 210433-3240 | Accolade II 132°, 40mm neck trial |
| <b>2</b>  | 210433-3237 | Accolade II 132°, 37mm neck trial |
| <b>3</b>  | 210433-3235 | Accolade II 132°, 35mm neck trial |
| <b>4</b>  | 210433-3230 | Accolade II 132°, 30mm neck trial |
| <b>5</b>  | 210433-3227 | Accolade II 132°, 27mm neck trial |
| <b>6</b>  | 210433-2737 | Accolade II 127°, 40mm neck trial |
| <b>7</b>  | 210433-2740 | Accolade II 127°, 37mm neck trial |
| <b>8</b>  | 210433-2735 | Accolade II 127°, 35mm neck trial |
| <b>9</b>  | 210433-2730 | Accolade II 127°, 30mm neck trial |
| <b>10</b> | 210433-2727 | Accolade II 127°, 27mm neck trial |
| <b>11</b> | 209870      | Anato neutral neck trial, left    |
| <b>12</b> | 209873      | Anato anteverted neck trial       |

**Description**

|           |                  |   |
|-----------|------------------|---|
| <b>13</b> | 209871           | Anato neutral neck trial, right                             |
| <b>14</b> | 111157           | Fixed screwdriver handle                                    |
| <b>15</b> | 110162           | Square driver   |
| <b>16</b> | 112730           | 3.5Mm hex driver  |
| <b>17</b> | 209830           | Mako inline-offset shell impactor (For trident shells only) |
| <b>18</b> | 212760           | Mako offset reamer handle (underneath impactor)             |
| <b>19</b> | 111120           | Bone pin adaptor  |
| <b>20</b> | 206270           | Mako shell impaction platform                               |
| <b>21</b> | 112340 or 110100 | Checkpoint driver   |
| <b>22</b> | 208074           | Mako straight reamer handle (underneath impactor)           |
| <b>23</b> | 209820           | Mako straight shell impactor                                |



- |          | Description    |
|----------|----------------|
| 1 209063 | MICS Handpiece |

---

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker's product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. The products depicted are CE marked according to the Medical Device Directive 93/42/EEC. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your sales representative if you have questions about the availability of Stryker's products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Accolade, Anato, Mako, Stryker, Trident. All other trademarks are trademarks of their respective owners or holders.